

Introducing The FORS™-15 Off-Loading Insole¹

The Durable and Adaptable Foot Off-Loading Solution

For Surgical Shoes

Closed-Toe



Darco® APB™



Darco® MedSurg™

Open-Toe



For Therapeutic Boots
(Boot not included)



4 Available Sizes, Trim-to-Fit Others

Consistent Off-Loading & Cushioning

- Durable memory resistant Poron® for consistent cushioning, crush resistance, and long life.
- Highly stable off-loading zones: Pre-cut plugs remain solidly in place unless removed. No unstable and shifting pegs means your patients will be consistently off-loaded.
- Effective offloading anywhere on the plantar surface, including along edges of the foot.
- Unique fabric mid-layer moderates “edge pressure” and prevents “bottoming out”.

Best in Class Materials

- Total of 15mm thick Poron® padding, with 5 mm always present under the foot. Multilayered design provides softness and support where needed: soft "slow rebound" Poron® adjacent to the foot, and denser Poron® below for support.
- Tough and abrasion resistant Alcantara® top cover minimizes slip and shear, absorbs perspiration, and is non-allergenic and bacteria resistant.

Uniquely adaptable:

- The FORS™-15 insole can be trimmed to any desired shape or size.

Easy Care:


- To completely clean and disinfect, soak the FORS™-15 insole in 5% - 10% chlorine bleach solution for 5 - 10 minutes and lightly scrub as needed, then thoroughly rinse and air dry.



www.bonapeda.com
(973) 415-9090
customercare@bonapeda.com

Bonapēda
PRODUCTS FOR HEALTHIER FEET

The FORS™-15 Off-Loading Insole¹ - Construction



FORS only (Open-Toe)	Shoe Size	Darco® MedSurg™ Size ²	Darco® Cat#	FORS with Darco Shoe
M494-15PZ (XL)	M 12.5 - 14	Men's X-Large	MQM4B	M495-SHOE-MX
M494-15PZ (LG)	M 10.5 - 12	Men's Large	MQM3B	M495-SHOE-ML
M494-15PZ (MED)	M 8.5 - 10	Men's Medium	MQM2B	M495-SHOE-MM
M494-15PZ (SM)	M 6-8	Men's Small	MQM1B	M495-SHOE-MS
M494-15PZ (W-L)	W 8.5 - 10	Women's Large	MQW3B	M495-SHOE-WL
M494-15PZ (W-M)	W 6.5 - 8	Women's Medium	MQW2B	M495-SHOE-WM
M494-15PZ (W-S)	W 4 - 6	Women's Small	MQW1B	M495-SHOE-WS

FORS Only (Closed-Toe)	Shoe Size	Darco® APB™ Size ²	Darco® Cat#	FORS + Darco Shoe
M494-15PZ-CT (XL)	W11.5 - 13 / M10.5 - 12.5	Extra-Large	APQ4B	M495-SHOE-CT-XL
M494-15PZ-CT (LG)	W 9.5 - 11 / M 8.5 - 10	Large	APQ3B	M495-SHOE-CT-L
M494-15PZ-CT (M)	W 8 - 9 / M 7 - 8	Medium	APQ2B	M495-SHOE-CT-M
M494-15PZ-CT (SM)	W 6 -7.5 / M 5 - 6.5	Small	APQ1B	M495-SHOE-CT-S
M494-15PZ-CT (XS)	W 3 - 5.5	Extra-Small	APQ0B	M495-SHOE-CT-XS

FORS Boot Insole	Shoe Size	DJO® CAM Boots ³	DJO® Size	(Boots Not Included)
M494-15PZ-B (SM)	< M4 / W5	DJO MaxTrax	SMALL	N/A
M494-15PZ-B (MED)	M4.5-6.5 / W5.5-10	DJO MaxTrax	MEDIUM	N/A
M494-15PZ-B (LG)	M6.5-11 / W10.5 - 13.5	DJO MaxTrax	LARGE	N/A
M494-15PZ-B (XLG)	M11.5-14.5 / W14 - 17	DJO MaxTrax	X-Large	N/A

1. U.S. Patent 10,441,031

2. Darco®, MedSurg™ and APB™ are trademarks of DARCO International, Inc., 810 Memorial Blvd, Huntington, WV 25701, USA

3. DJO®, DONJOY® and MaxTrax™ are trademarks of DJO, Global, 1430 Decision Street, Vista, CA 92081, USA



www.bonapeda.com
 (973) 415-9090
customercare@bonapeda.com

Bonapēda
 PRODUCTS FOR HEALTHIER FEET